

### WHAT ARE MEDICAL COMPRESSION GARMENTS (STOCKINGS)

Compression stockings are specialized garments made with special fabric. The importance of compression socks lies in their ability to apply gentle pressure to the legs, enhancing venous return. Primarily Compression Stockings have two functions: -

• PREVENTION OF VENOUS DISORDERS VISTA CARE Comfort That Heals

• MANAGEMENT OF VENOUS DISORDERS

This is achieved mainly by providing desired compression thereby helping in:

- IMPROVED BLOOD CIRCULATION
- REDUCE SWELLLING (EDEMAS)









## WHO SHOULD USE MEDICAL STOCKINGS (PREVENTIVE)

- <u>ELDERLY POPULATION</u> for chronic venous insufficiency, varicose veins & DVT prevention
- PREGNANT WOMEN Relief from Varicose Veins and swelling




## WHO SHOULD USE MEDICAL STOCKINGS (PREVENTIVE)

- LONG WORK-HOURS PROFESSIONALS Prevents Varicose Veins & Swelling,
   reduce muscle soreness
- ATHLETES Prevention of Varicose Veins, swelling and enhanced performance due to confort












## WHO SHOULD USE MEDICAL STOCKINGS (MANAGEMENT)

Vista Care Medical Compression Garments are useful in the treatment of:

- Varicose Veins
- Venous Ulcers
- Venous insufficiency
- Lymphoedema
- Deep vein Thrombosis (DVT)
- Cellulitis
- Keloids
- Postural Hypotension
- Post Liposuction
- Post Mastectomy Lymphoedema
- Abdominal Surgeries





DVT



### STOCKINGS (MANAGEMENT)

5

**VARICOSE VEINS** 



- **POST SURGICAL CARE** Anti embolism therapy and post operative recovery
- <u>VENOUS DISORDERS</u> Improved blood flow, reduce swelling



**BURN KELOIDS** 



**VENOUS ULCERS** 



WHO SHOULD USE MEDICAL

LYMPHOEDEMA




### WHY VISTA CARE?


Specially developed Fabric that is: Comfort That Heals


- All Direction Stretchable Thus providing comfort fit.
- Cotton-Rich, Cool Fabric that breathes
- Water repellant, anti-microbial & Skin safe
- Especially designed for Vietnam Climate and mobility needs.

Graduated Compression for better blood flow

Medical Compression Stockings (Class I & II)








## VISTA CARE COMPRESSION STOCKINGS - MODELS:

Vista care Medical Stockings are available in three major models:

BELOW KNEE

• THIGH HIGH

• FULL LEG




# VISTA CARE COMPRESSION STOCKINGS - MODELS & SIZE CHART:

Vista care Medical Stockings are available in three major models:



#### STANDARD SIZE MEASUREMENTS



All measurements (cm)		XS	S	M	L	XL
KNEE LENGTH	Min Ankle (c4)	16-18	18-21	21-24	24-28	28-32
	Max Calf (c5)	24-28	28-34	34-38	38-42	42-46
	Just below knee (c6)	22-24	26-28	28-30	30-32	32-35
	Length (F)	28-30	30-32	32-34	34-36	36-39
	Min Ankle (c4)	16-18	18-21	21-24	24-28	28-32
THIGH	Max Calf (c5)	24-28	28-34	34-38	38-42	42-46
HIGH	Just below knee (c6)	22-24	26-28	28-30	30-32	32-35
LENGTH	Mid of Thigh (c9)	36-40	40-44	44-48	48-52	52-56
	Length (C+D)	56-60	60-64	64-68	68-72	72-76
	Min Ankle (c4)	16-18	18-21	21-24	24-28	28-32
	Max Calf (c5)	24-28	28-34	34-38	38-42	42-46
FULL LEG	Just below knee (c6)	22-24	26-28	28-30	30-32	32-35
LENGTH	Mid of Thigh (c9)	36-40	40-44	44-48	48-52	52-56
	Gluteal Fold (c10)	44-50	50-54	54-58	58-62	62-66
	Length (B+D)	60-64	64-68	68-72	72-76	76-80



# WHY VISTA CARE COMPRESSION STOCKINGS (MANAGEMENT):

- <u>Custom-size</u> to exact size-in the condition of the disease
- A limb with inflammation / Edema / swelling needs <u>appropriate</u>
   <u>compression</u> at all the points
- Hence, Vista care garment is prepared to exact size at any stage of the disease.
- Progressive correction is part of the after sales service when swelling recedes








## VISTA CARE COMPRESSION GARMENTS - ARMS & BODY SLEEVES:

- Vista care compression garments
  have wide applications for arms
  and other parts of the human
  body, as it is custom -made and
  Engineered to provide
  compression & support.
- Also for the treatment of Post-Burn Hypertrophies.

### **VISTAAR COMPANY POLICY:**

- Individual attention & personal care
- Trained staff visits patient to take the <u>measurements at all the requisite</u>
   <u>points</u> with utmost care
- Garment manufactured to provide the <u>appropriate graduated compression</u> at various points based on patient's condition
- Garment prepared within 1-2 working days
- Progressive correction provided (part of after sales service)



VISTAAR INDIA CG CO.
LTD.,



VISTAAR INDIA CG CO.
LTD.,

### **GET IN TOUCH**

- Mumbai, Maharashtra (INDIA)
- +91 9702274489
- +91 9702581163
- New Delhi, INDIA
- +91 9958448607 / +91 9818291300
- vistaarinternationalcoltd@gmail.com



VISTAAR INDIA CG CO.
LTD.,

### COMING SOON...

- O HO CHI MINH CITY,
  VIET NAM
- +84-909774302
- +91 9702404473
- vistaarinternationalcoltd@gmail.com

https://www.vistacareindia.com