



VISTA CARE  
Comfort That Heals

# MEDICAL COMPRESSION GARMENTS



VISTAAR INTERNATIONAL CG CO. LTD.,

<https://www.vistacareindia.com>

# MEDICAL COMPRESSION STOCKINGS

**VISTA CARE** provides advanced **Medical Compression Therapy**, the definitive non-invasive standard for the management of venous and lymphatic disorders.

## **THERAPEUTIC MECHANISM & OBJECTIVES:**

**Precision Graduation:** Specifically engineered fabric that apply controlled, graduated pressure to restore healthy hemodynamics.

**Dual-Action** Therapy: Validated for both the prophylactic prevention and clinical management of vascular pathologies.

**Optimized Venous Return:** Enhances blood flow velocity and prevents venous stasis.

**Edema Mitigation:** Effectively reduces interstitial fluid accumulation and manages limb swelling.

**VISTA CARE** manufactures medical-grade compression garments compliant with international standards: **EN 13770, DIN 58133, and RAL-GZ 387**. These standards ensure that every garment delivers precise, consistent pressure profiles optimized for clinical efficacy.






# FEATURES OF VISTA CARE STOCKINGS

- **GRADUATED COMPRESSION**


VISTA CARE medical compression stockings are engineered to provide graduated compression to the affected limb so as to enable effective blood flow.



- **COMFORT THAT HEALS**

**VISTA CARE** medical compression stockings are designed with specially developed cotton fabric that is

- > **All Directional Stretch**
- > **Anti-microbial & non-irritating**
- > **Breathable & Soft to Skin**
- > **Long Durability**
- > **Comfortable to wear & maintain in all weather conditions**




# WHY VISTA CARE MEDICAL COMPRESSION STOCKINGS ?

Feature	Medical-Grade	OTC Products
<b>Pressure Accuracy</b>	±5% tolerance per RAL-GZ 387	No standard testing
<b>Pressure Profile</b>	Graduated compression (precise mmHg values)	Variable, inconsistent
<b>Clinical Testing</b>	Tested with pressure sensors (HOSYcan)	No formal testing
<b>Standards Compliance</b>	EN 13770, DIN 58133, RAL-GZ 387	None
<b>Clinical Efficacy</b>	Proven in clinical trials	Not validated
<b>Durability</b>	6+ months of therapeutic use	3-4 months typical
<b>Customisation</b>	Individual fit and pressure profiles	Limited options
<b>Medical Device Status</b>	Regulated medical device	Consumer textile



# WHO SHOULD USE VISTA CARE Class I & II (18-32mmHg) STOCKINGS (PREVENTIVE)

**Elderly Population** for chronic venous insufficiency, Varicose Veins & DVT



**Pregnant women** throughout pregnancy for prevention / treatment of Varicose Veins and Edemas

**Long work-hour professionals -**  
Prevention of Varicose Veins,  
Swelling and muscle soreness




**Athletes** for prevention of Varicose Veins, Edemas and enhanced performance due to muscle comfort



# STANDARD SIZE MEASUREMENTS CHART FOR APPROPRIATE SIZE SELECTION

STANDARD SIZE MEASUREMENTS						
All measurements (cm)		XS	S	M	L	XL
<b>BELOW KNEE</b>	Min Ankle (cB)	16-19	20-22	23-25	26-28	29-31
	Max Calf (cC)	27-36	30-39	33-41	36-44	39-47
	Just below	26-35	29-38	32=40	35-43	38-46
	Length (LD)	30-32	32-34	34-36	36-39	36-39
<b>THIGH HIGH</b>	Min Ankle (cB)	16-19	20-22	23-25	26-28	29-31
	Max Calf (cC)	27-36	30-39	33-41	36-44	39-47
	Just below	26-35	29-38	32=40	35-43	38-46
	Mid of Thigh (cF)	40-44	44-48	48-52	52-56	56-62
	Length (LF)	58-62	62-64	62-64	62-64	62-64
<b>FULL LEG</b>	Min Ankle (cB)	16-19	20-22	23-25	26-28	29-31
	Max Calf (cC)	27-36	30-39	33-41	36-44	39-47
	Just below	26-35	29-38	32=40	35-43	38-46
	Mid of Thigh (cF)	40-44	44-48	48-52	52-56	56-62
	Gluteal Fold (cG)	44-59	49-63	54-66	58-70	63-74
	Length (LG)	64-68	68-72	68-72	68-72	68-72



# MANAGEMENT OF VENOUS / LYMPHATIC DISORDERS (CUSTOMIZED SOLUTIONS)

VISTA CARE medical stockings  
MANUFACTURES  
**CUSTOMIZED STOCKINGS**  
THEREBY PROVIDING  
EFFECTIVE SOLUTION BY:

- TAKING **EXACT MEASUREMENTS** OF THE PATIENT / CONDITION **AT THE DOOR STEP**
- PREPARING THE REQUISITE GARMENT & DELIVERY WITHIN SHORT TURN-AROUND TIME
- **PROGRESSIVE CORRECTION** PROVIDED FOR PATIENTS REQUIRING RE-SIZING OF THE GARMENT POST RECOVERY



Chronic Vein Insufficiency



Varicose Veins



Lymphoedema



Venous Ulcers



Post Thrombotic Syndrome



Edema



Dee Vein Thrombosis (DVT)




Cellulitis



Pregnancy related issues

# VISTA CARE PRODUCT RANGE



**Below Knee Stockings**



**Thigh High Stockings**



**Full Pantyhose**



**Full Leg Stockings**



**Armsleeve**



**Chin Support**



**Full Vest Garment**




**Torso**



**Bottom Girdle**



**Full Body Suit**



# GET IN TOUCH



**Mumbai, Maharashtra (INDIA)**



+91 9702274489



+91 9702581163



**New Delhi, INDIA**



+91 9958448607 / +91 9818291300



**Ho Chi Minh City, VIET NAM**



+91 9702404473



+84 909774302



vistaarinternationalcoltd@gmail.com

<https://www.vistacreindia.com>

